


# **TECHNICAL DATA SHEET**



# **PRODUCT SPECIFICATION**

Flow movement	Hybrid axial/radial
Flow direction	Down
Weight	8.00 ± 0.25 g
Viscosity range	1-205 cP
RPM range	100-1100 RPM
Volume range	100-500 mL
Dimension	40 (L) x 13 (W) x 15 (H) mm

# SAN-01



#### Preferred vessel:



**ROTATIONAL SPEED (RPM)** 

# **HIGH SHEAR DESIGN**



#### Utilized in emulsion and dispersion processes of :

- Emulsion and suspension polymerization
- Cosmetic creams, ceramic slurries, and polishing compounds
- Conductive inks and nanoparticle suspensions

## **MATERIAL TECHNICAL DATA**


Shell material	polyethylene naphthalate (PEN)  🎪
Application temperature	-60 to 155°C
Color	white
Coefficient of friction	0.12
Hydrolysis resistance	200 hours
Compliance	FDA food contact, USP <87>
Magnet material	SmCo permanent
Remanence (Br)	9,800~11,000 Gs
Sterilization	gamma, steam, EtOH

# PRODUCT HIGHLIGHT

PFAS FREE	Shell made of PEN to reduce our dependency on PTFE & PFAS
	Food contact safe
USD. < 87 >	No cell toxicity, bio-inert
ROHS2 2011/65/EU	Non-hazardous materials
5°CO2	Estimated carbon footprint of PEN is $3.4 \text{ kg CO}_2 \text{ eq/kg}$
	Recyclable with PET bottle

## **FLOW DIRECTION SETTING**


The up/down flow direction depends on the spin direction of the plate and the handedness of the product.



## RECOMMENDED OPERATING RANGE

.

K




# **CHEMICAL RESISTANCE CHART**

	PEN	PVDF	PTFE	
Acetic acid (10%)	А	А	А	
Acetic acid (40%)	А	А	А	
Acetone	В	D	А	
Ammonium chloride (25%)	Α	Α	А	
Benzene	В	Α	А	
Chloroform	С	Α	А	
Cyclohexane	Α	Α	А	
Cyclohexanol	Α	Α	А	
Calcium chloride (10%)	Α	Α	А	
Chromic acid (10%)	Α	Α	А	
Citric acid	Α	Α	А	
Diethyl ether	Α	Α	А	
Dimethyl sulfoxide (DMSO)	С	С	А	
Dimethyl formamide (DMF)	С	D	А	
Ethanol	Α	Α	А	
Ethyl acetate	Α	D	А	
Ethylene diamine	С	В	А	
Ethylene glycol	Α	Α	А	
Formaldehyde (40%)	Α	Α	А	
Formic acid (10%)	Α	Α	А	
Formic acid (30%)	С	В	А	
Gasoline	А	А	А	
Glycerol	А	А	А	
Heptane	А	Α	А	
Hexane	Α	Α	А	
Hydrobromic acid (10%)	А	А	А	
Hydrobromic acid (47%)	А	А	А	
Hydrochloric acid (10%)	А	А	А	
Hydrochloric acid (37%)	А	А	А	
Hydrofluoric acid (5%)	А	А	А	
Hydrofluoric acid (50%)	D	В	А	

	PEN	PVDF	PTFE
Hydrogen peroxide (30%)	А	А	А
Isopropyl alcohol (IPA)	А	А	А
Methanol	А	А	А
Methyl ethyl ketone (MEK)	А	D	А
Nitric acid (10%)	А	А	А
Nitric acid (20%)	А	А	А
Nitric acid (67%)	D	А	А
Phosporic acid (10%)	А	В	А
Phosphoric acid (85%)	D	В	А
Potassium hydroxide (10%)	А	А	А
Potassium hydroxide (30%)	В	А	А
Sulfuric acid (10%)	А	А	А
Sulfuric acid (30%)	С	А	А
Sulfuric acid (>80%)	D	А	А
Sodium acetate (40%) sol.	А	А	А
Sodium chloride (sat.) (32%)	А	А	А
Sodium bicarbonate (sat.)	А	А	А
Sodium hydroxide (10%) sol.	А	А	А
Sodium hydroxide (30%) so <b>l</b> .	А	С	А
Sodium hypochlorite	А	А	А
Styrene (monomer)	А	А	А
Tetrahydrofuran (THF)	В	В	А
Toluene	А	А	А
Triethanolamine	С	А	А
Vinyl chloride (monomer)	А	А	А
Xylene	А	А	А
Resistance at 23°C, immersed for 30 da Ranking definition: A: weight change <1%, tensile > 95% B: weight change between 1~10%, ter C: weight change >10%, tensile <75% D: dissolved or swelled		ó	

## **SELECTION MATRIX**

type of magnetic stir bars categorized by performance toward various speeds and viscosities.



## PRECAUTION

#### Hazard Identification:

Not classified as hazardous according to OSHA or GHS regulations. Avoid exposure to DMSO, DMF, CHCl<sub>a</sub>, CCl<sub>a</sub>, and HF.

#### Check compatibility:

Use chemical resistance chart as a reference.

### Conditions to Avoid:

Temperatures above 160°C.

#### Storage:

Magnetic stir bars should be stored in pairs to maintain their magnetic strength, away from steel or iron surfaces and other magnetic fields. Keep them at ambient temperatures to avoid thermal stress, and handle carefully to prevent mechanical damage.

#### **Disposal Considerations:**

PEN (Polyethylene Naphthalate) can be recycled alongside PET (Polyethylene Terephthalate). Recycled as much as possible with PET or through incineration.

## FLXR Engineering Co., Ltd.

♀ 1F & 2F, 11 - 1 Wuquan 1<sup>#</sup> Rd., Xinzhuang Dist., New Taipei City, Taiwan 24892
▶ +886-2-2290-1122
₱ +886-2-2299-5222
⊕ hello@flxr.engineering





